

Parkinson's rollator with
customizable cues

rollz motion
rhythm

Regain your walking rhythm

- Laser, sound and vibration cues
- Great stability and comfort
- Advanced and smart cue control
- Easy to push and steer
- Supports various rehabilitation programs

Designed in the Netherlands



Step out of a freeze with the Rollz Motion Rhythm

This stable rollator is developed for people with a disrupted gait pattern or who suffer from freezing episodes. It encourages the start of a rhythmic gait pattern with the help of three cues.

Innovation

Specialized version of the Rollz Motion that addresses people with Parkinson’s or a disrupted gait pattern.

Three cues

Laser line on the ground, metronome sound and vibration in the handlebars.



Laser



Vibration



Sound

Adjustable

Settings can be customized in accordance with the rehabilitation program.

Start the rhythm

The desired cues combination and tempo are easy to set and activate.

Stable

Offers unmatched stability to prevent falls that could occur due to balance problems.

Smartphone app


Allows to create random patterns for preventing habituation.



Accessorize the Parkinson’s rollator



Wheelchair package



Drag brakes



Cup holder



Light



Back support



Seat belt



Travel cover



Locking hook

Frame colour



Matt black

Specifications

Rollator weight	11,6kg
Maximum weight capacity	125kg
User height	1,60 - 1,90m
Seat height	55cm
Handle height range	87 - 98cm
Overall width x length	67 x 63cm
Width folded	31cm